

[x](#)

How to Use

- **Searching:** Type keyword in search field at top of page. Search by all or part of a monograph title. For searches using multiple criteria, you will find items that match each of the specified criteria unless quotation marks are used.
 - For example, a search on Aminosalicyclic Acid Tablets will result in anything that contains “Aminosalicyclic” OR “Acid” OR “Tablets”
 - A search for “Aminosalicyclic Acid Tablets” will result in anything that specifically contains “Aminosalicyclic Acid Tablets”
- **Sorting:** Click on any column header title to sort alphabetically or chronologically in ascending or descending order. Note: the page load column is sorted alphabetically so that a number is ordered by first digit vs. by the actual number; thus, numbers will not always be in order.
 - For example, page 2178 will come before page 74 on a page sort.
- **Downloading:** You can download the Errata table in Comma-separated Value (.csv). The download will include the Errata that you have filtered on.
- **Importing:** You will need to import the file into Excel or Open Office with UTF-8 encoding, as opposed to simply opening it. To import, open Excel or Open Office and select import from the File drop-down. Depending on the version you are using, you should be presented with import formatting options to include UTF-8 as one of the first steps. Importing via UTF-8 should eliminate odd character conversions.

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
LEVOFLOXACI ADDITIONAL R	USP36–NF31	4099	26-Jul-2013	1-Aug-2013	USP38–NF33	First	Line 2 of USP

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
N	EQUIREMENT S/USP Reference Standards <11>							Supplement to USP37–NF32	Levofloxacin Related Compound A RS: Change (S)-9-Fluoro-3-me thyl-10-(piperazi n-1-yl)-7-oxo-2, 3-di hydro-7H -p yrido [1,2,3- <i>de</i>][1,4-benzoxazi ne-6-carbocyclic acid. to: (S)-9-Fluoro-3-me thyl-10-(piperazi n-1-yl)-7-oxo-2, 3-di hydro-7H -p yrido [1,2,3- <i>de</i>][1,4]benzoxazi ne-6-carboxylic acid. AND Line 2 of USP

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
VENLAFAXINE ADDITIONAL R HYDROCHLOR EQUIREMENT	USP36–NF31	5551	26-Jul-2013	1-Aug-2013	USP38–NF33	First Supplement to	Levofloxacin Related Compound B RS: Change (S)-9,10-Difluoro-3-methyl-7-oxo-2,3-di hydro-7H -p yrido [1,2,3- <i>de</i>][1,4-benzoxazi ne-6-carbocylic acid. to: (S)-9,10-Difluoro-3-methyl-7-oxo-2,3-di hydro-7H -p yrido [1,2,3- <i>de</i>][1,4]benzoxazi ne-6-carboxylic acid. Line 2 of Venlafaxine

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Errata Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
IDE	<i>SIUSP Reference Standards <11></i>							<i>USP37–NF32</i>	Related Compound A RS: Change 1-(1-(4-methoxyphenyl)-2-(meth ylamino)ethyl)cy clohexanol. $C_{16}H_{25}NO_2$ 263.38 to: 1-(1-(4-Methox yphenyl)-2-(met hylamino)ethyl) cyclohexanol hydrochloride. $C_{16}H_{25}NO_2 \cdot$ HCl 299.84
STINGING NETTLE	COMPOSITION <i>/Content of Total Amino Acids</i>	<i>USP36–NF31</i>	1604	26-Jul-2013		1-Aug-2013	<i>USP38–NF33</i>	<i>First Supplement to USP37–NF32</i>	Line 1 of <i>Reagent solution:</i> Change Solution containing 1.00 g of ninhydrin, 1.50 g of hydrindantin, to: Solution containing 1.00 g of ninhydrin,

Monograph Title Section		Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
GENTAMICIN SULFATE	IM PURITIES/ <i>Limit of Methanol</i>	<i>First Supplement to USP36–NF31</i>	5990	26-Jul-2013	1-Aug-2013	USP38–NF33	<i>First Supplement to USP37–NF32</i>	150 mg of hydrindantin, Line 11 of <i>Analysis:</i> Change Result = $(R_U/R_S) \times (C_S/C_U) \times D \times F \times 100$ to: Result = $(R_U/R_S) \times (C_S/C_U) \times D \times F$ AND Line 18 of <i>Analysis:</i> Change C_S = percentage of methanol in the <i>Standard solution</i> (% v/v) to: C_S = percentage of methanol in the <i>Standard solution</i> AND Line 23 of <i>Analysis:</i>

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Errata Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
AMOXICILLIN	IM PURITIES/ <i>Organic Impurities/Procedure</i>	USP36–NF31	2477	26-Jul-2013		1-Aug-2013	USP38–NF33	<i>First Supplement to USP37–NF32</i>	Change <i>F</i> = conversion factor, 0.001 g/mg to: <i>F</i> = conversion factor, 1000 mg/g Line 2 of <i>Acceptance criteria</i> : Change [Note—The reporting limit is 0.03% of the amoxicillin peak from the <i>Standard solution</i> .] to: [Note—The reporting limit is 0.03 times the amoxicillin peak from the <i>Standard solution</i> .]
CLARITHROMYCIN EXTENDED-RELEASE TABLETS	PERFORMANCE TESTS/ <i>Dissolution <711>/Test 4</i>	USP36–NF31	3019	26-Jul-2013		1-Aug-2013	USP38–NF33	<i>First Supplement to USP37–NF32</i>	Line 2 of <i>Standard solution</i> : Change and Medium

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
DIPHENHYDRAMINE HYDROCHLORIDE ORAL SOLUTION	USP36–NF31	3277	26-Jul-2013	1-Aug-2013	USP38–NF33	First Supplement to USP37–NF32	(96:4). to: and <i>Medium</i> (4:96). Line 2: Change <i>Mobile phase, Standard preparation, System suitability solution, and Chromatographic system</i> —Prepare as directed in the Assay under <i>Diphenhydramine Hydrochloride</i> . to: <i>Mobile phase</i> —Prepare a solution of acetonitrile, water, and triethylamine (50: 50: 0.5), adjust with glacial acetic acid to a pH of 6.5, filter, and

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
							degas. Make adjustments if necessary (see <i>System Suitability</i> under <i>Chromatography</i> <621>). <i>Standard preparation</i> —Dissolve an accurately weighed quantity of USP Diphenhydramine Hydrochloride RS in water to obtain a solution having a known concentration of about 0.5 mg per mL. AND After the <i>Assay preparation</i> subsection: Add <i>System suitability solution</i>

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
							—Dissolve about 5 mg of benzophenone in 5 mL of acetonitrile, dilute with water to 100 mL, and mix. Transfer 1.0 mL of this solution and 5 mg of diphenhydramine hydrochloride to a 10-mL volumetric flask, dilute with water to volume, and mix. <i>Chromatographic system (see Chromatography <621>)</i>—The liquid chromatograph is equipped with a 254-nm detector and a 4.6-mm × 25-cm column that contains packing L10.

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
							The flow rate is about 1 mL per minute. Chromatograph the <i>System suitability solution</i>, and record the peak responses as directed for <i>Procedure</i>; the resolution, <i>R</i>, between the benzophenone and diphenhydramine peaks is not less than 2.0. Chromatograph replicate injections of the <i>Standard preparation</i>, and record the peak responses as directed for <i>Procedure</i>; the relative standard deviation is not more than

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
							2.0%, and the tailing factor for the diphenhydramine hydrochloride peak is not more than 2.0. AND Line 1 of <i>Procedure:</i> Change Proceed as directed for <i>Procedure</i> in the Assay under <i>Diphenhydramine Hydrochloride</i>. to: Separately inject equal volumes, about 10 µL, of the <i>Standard preparation</i> and the Assay <i>preparation</i> into the chromatograph, record the chromatograms,

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
FLUPHENAZIN Assay E DECANOATE INJECTION	USP36–NF31	3639	26-Jul-2013	1-Aug-2013	USP38–NF33	First Supplement to USP37–NF32	and measure the responses for the major peaks. Line 11 of <i>Standard preparation</i> : Delete (1:5) AND Line 8 of <i>Assay preparation</i> : Delete (1:5)
OLMESARTAN IM MEDOXOMIL PUR ITIES/ <i>Organic Impurities/Impurity Table</i>	USP36–NF31	4570	26-Jul-2013	1-Aug-2013	USP38–NF33	First Supplement to USP37–NF32	Footnote d: Change ((5-Methyl-2-oxo-1,3-dioxol-4-yl)methyl 4-(2-hydroxypropan-2-yl)-2-propyl-1-((2'-(1-trityl-1 <i>H</i> -tetrazol-5-yl)biphenyl-4-yl)methyl)-1 <i>H</i> -imidazole-5-carboxylate. to: (5-Methyl-2-oxo-1,3-dioxol-4-yl

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Errata Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
									)methyl 4-(2-hydroxypropan-2-yl)-2-propyl-1-((2'-(2-trityl-1H-tetrazol-5-yl)biphenyl-4-yl)methyl)-1H-imidazole-5-carboxylate.
POWDERED GYMNEMA	IDENTIFICATION/ N/B. Thin-Layer Chromatography/ Chromatographic system	<i>First Supplement to USP36–NF31</i>	5883	26-Jul-2013		1-Aug-2013	<i>USP38–NF33</i>	<i>First Supplement to USP37–NF32</i>	Line 2 of Adsorbent: Change 5m to: 5 µm
GLYCERYL DISTEARATE	ASSAY/ Procedure/ Chromatographic system	<i>USP36–NF31</i>	2029	26-Jul-2013		1-Aug-2013	<i>USP38–NF33</i>	<i>First Supplement to USP37–NF32</i>	Line 1 of Column temperature: Change Column temperature: 40° to: Temperatures Detector: 40° Column: 40°
OXCARBAZEPINE	IDENTIFICATION/ Purities/Organic Impurities,	<i>First Supplement to USP36–NF31</i>	6035	26-Jul-2013		1-Aug-2013	<i>USP38–NF33</i>	<i>First Supplement to USP37–NF32</i>	Row 9 of Column 1 of Table 3: Change

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
	Procedure 2								Oxcarbazepine related compound E ^g to: Oxcarbazepine related compound E AND Delete footnote g
CEFDINIR CAPSULES	IMPURITIES/Organic Impurities/Table 2	USP36–NF31	2850	26-Jul-2013		1-Aug-2013	USP38–NF33	First Supplement to USP37–NF32	Row 16 of Column 1: Change Cefdinir impurity 2 ^e to: Cefdinir impurity 2 ^f AND Row 21 of Column 1: Change Cefdinir impurity 3 ^e to: Cefdinir impurity 3 ^f
DIPHENHYDRAmine HYDROCHLORIDE CAPSULES	Assay	USP36–NF31	3276	26-Jul-2013		1-Aug-2013	USP38–NF33	First Supplement to USP37–NF32	Line 2: Change Mobile phase, Standard preparation,

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
							<i>System suitability solution</i>, and <i>Chromatography</i> c <i>system</i>—Prepare as directed in the Assay under <i>Diphenhydramine Hydrochloride</i>. to: <i>Mobile phase</i>—Prepare a solution of acetonitrile, water, and triethylamine (50: 50: 0.5), adjust with glacial acetic acid to a pH of 6.5, filter, and degas. Make adjustments if necessary (see <i>System Suitability</i> under <i>Chromatography</i> <621>). <i>Standard</i>

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
							<i>preparation</i> —Dissolve an accurately weighed quantity of USP Diphenhydramine Hydrochloride RS in water to obtain a solution having a known concentration of about 0.5 mg per mL. AND After the Assay <i>preparation</i> subsection: Add <i>System suitability solution</i> —Dissolve about 5 mg of benzophenone in 5 mL of acetonitrile, dilute with water to 100 mL, and mix. Transfer

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
							1.0 mL of this solution and 5 mg of diphenhydramine hydrochloride to a 10-mL volumetric flask, dilute with water to volume, and mix. <i>Chromatographic system (see Chromatography <621>)</i>—The liquid chromatograph is equipped with a 254-nm detector and a 4.6-mm × 25-cm column that contains packing L10. The flow rate is about 1 mL per minute. Chromatograph the System suitability solution, and record the peak

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
							responses as directed for <i>Procedure</i>; the resolution, <i>R</i>, between the benzophenone and diphenhydr amine peaks is not less than 2.0. Chromatograph replicate injections of the <i>Standard preparation</i>, and record the peak responses as directed for <i>Procedure</i>; the relative standard deviation is not more than 2.0%, and the tailing factor for the diphenhydra mine hydrochloride peak is not more than 2.0. AND

Monograph Title Section		Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
								Line 1 of <i>Procedure:</i> Change Proceed as directed for <i>Procedure</i> in the Assay under <i>Diphenhydramine Hydrochloride</i> . to: Separately inject equal volumes, about 10 µL, of the <i>Standard preparation</i> and the Assay <i>preparation</i> into the chromatograph, record the chromatograms, and measure the responses for the major peaks.
FELBAMATE TABLETS	IMPURITIES/Organic Impurities	USP36–NF31	3537	26-Jul-2013	1-Aug-2013	USP38–NF33	First Supplement to USP37–NF32	Line 1 of <i>Resolution:</i> Change NMT 2

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
	<i>ties/System suitability/Suitability requirements</i>								to: NLT 2
MOXIFLOXACIN OPHTHALMIC SOLUTION	<i>Related compounds</i>	<i>USP36–NF31</i>	4414	26-Jul-2013		1-Aug-2013	<i>USP38–NF33</i>	<i>First Supplement to USP37–NF32</i>	Row 1 of Column 3 of <i>Table 2</i> : Change Relative Retention Time vs. Moxifloxacin to: Relative Retention vs. Moxifloxacin AND Row 1 of Column 3 of <i>Table 3</i> : Change Relative Retention Time vs. Moxifloxacin to: Relative Retention vs. Moxifloxacin
VENLAFAXINE TABLETS	ADDITIONAL REQUIREMENT <i>SI/USP Reference</i>	<i>USP36–NF31</i>	5554	26-Jul-2013		1-Aug-2013	<i>USP38–NF33</i>	<i>First Supplement to USP37–NF32</i>	Line 2 of Venlafaxine Related Compound A

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
<i>Standards <11></i>									RS: Change 1-(1-(4-methoxy phenyl)-2-(meth ylamino)ethyl)cy clohexanol. $C_{16}H_{25}NO_2$ 263.38 to: 1-(1-(4-Methox yphenyl)-2-(met hylamino)ethyl) cyclohexanol hydrochloride. $C_{16}H_{25}NO_2 \cdot$ HCl 299.84
POWDERED STINGING NETTLE	COMPOSITION <i>/Content of Total Amino Acids</i>	USP36–NF31	1606	26-Jul-2013		1-Aug-2013	USP38–NF33	<i>First Supplement to USP37–NF32</i>	Line 1 of <i>Reagent solution:</i> Change Solution containing 1.00 g of ninhydrin, 1.50 g of hydrindantin, to: Solution containing 1.00 g of ninhydrin, 150 mg of hydrindantin,

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
MEPROBAMATE TABLETS	ASSAY/ Procedure	<i>First Supplement to USP36–NF31</i>	6015	26-Jul-2013		1-Aug-2013	USP38–NF33	<i>First Supplement to USP37–NF32</i>	Line 4 of <i>Standard solution</i> : Change Dissolve in 30% of the final flask volume, and dilute with water to volume. to: Dissolve in 30% of the final flask volume of acetonitrile, and dilute with water to volume.
BETAMETHASONE SODIUM PHOSPHATE	<i>Identification/B. Thin-Layer Chromatographic Identification Test <201></i>	USP36–NF31	2645	26-Jul-2013		1-Aug-2013	USP38–NF33	<i>First Supplement to USP37–NF32</i>	Line 1 of <i>Test solution</i> : Change 1 mg per mL. to: 1 mg per mL in methanol.
DACARBAZINE FOR INJECTION	<i>USP Reference standards <11></i>	USP36–NF31	3137	26-Jul-2013		1-Aug-2013	USP38–NF33	<i>First Supplement to USP37–NF32</i>	Line 3 of USP Dacarbazine Related Compound B RS: Change $C_4H_3N_5O$ 137.10 to:

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Errata Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
EDETATE DISODIUM	ASSAY/ <i>Procedure</i>	USP36–NF31	3370	26-Jul-2013		1-Aug-2013	USP38–NF33	<i>First Supplement to USP37–NF32</i>	$C_4H_3N_5O \cdot H_2O$ 155.12 Line 5: Delete <i>Titrimetric system</i> (See <i>Titrimetry</i> <541>.) <i>Mode</i> : Direct titration <i>Titrant</i> : 0.1 N sodium hydroxide VS <i>Endpoint detection</i> : Visual
LAMOTRIGINE TABLETS	PERFORMANCE TESTS/ <i>Dissolution</i> <711>/ <i>Test 1</i>	USP36–NF31	4056	26-Jul-2013		1-Aug-2013	USP38–NF33	<i>First Supplement to USP37–NF32</i>	Line 3 of <i>Standard solution</i> : Change 0.028 µg/mL to: 0.028 mg/mL Footnote c: Change 1-(3-Methoxyphenyl)-2-(dimethylaminomethyl)cyclohex-1-ene hydrochloride (identified and
TRAMADOL HYDROCHLORIDE EXTENDED-RELEASE TABLETS	HYMPURITIES/Organic Impurities/ <i>Table 2</i>	USP36–NF31	5438	26-Jul-2013		1-Aug-2013	USP38–NF33	<i>First Supplement to USP37–NF32</i>	Footnote c: Change 1-(3-Methoxyphenyl)-2-(dimethylaminomethyl)cyclohex-1-ene hydrochloride (identified and

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
							reported as an individual unspecified impurity if present). to: 1-(3-Methoxyphenyl)-2-(dimethylaminomethyl)cyclohex-6-ene hydrochloride (identified and reported as an individual unspecified impurity if present). AND Footnote d: Change 1-(3-Methoxyphenyl)-2-(dimethylaminomethyl)cyclohex-6-ene hydrochloride (identified and reported as an individual unspecified impurity if present).

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Errata Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
PURIFIED GYMNEMA EXTRACT	COMPOSITION <i>/Content of Gymnemic Acids</i>	<i>First Supplement to USP36–NF31</i>	5884	26-Jul-2013		1-Aug-2013	<i>USP38–NF33</i>	<i>First Supplement to USP37–NF32</i>	to: 1-(3-Methoxyphenyl)-2-(dimethylaminomethyl)cyclohex-1-ene hydrochloride (identified and reported as an individual unspecified impurity if present). Line 1 of <i>Acceptance criteria</i> : Change 90%–110% of the labeled amount to: 90.0%–110.0% of the labeled amount on the dried basis
SODIUM HYDROXIDE	ASSAY/ <i>Procedure</i>	<i>USP36–NF31</i>	2203	26-Jul-2013		1-Aug-2013	<i>USP38–NF33</i>	<i>First Supplement to USP37–NF32</i>	Line 11 of <i>Analysis</i> : Change = volume of <i>Titrant</i> consumed by the <i>Sample</i> to

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
CLARITHROM YCIN FOR ORAL SUSPENSION	ASSAY/ <i>Procedure</i>	USP36–NF31	3018	26-Jul-2013		1-Aug-2013	USP38–NF33	<i>First Supplement to USP37–NF32</i>	the first endpoint (mL) to: = volume of <i>Titrant</i> consumed by the <i>Sample</i> to the second endpoint (mL) Change the subsection head <i>Buffer</i>: to: <i>Buffer A</i>: AND After the <i>Buffer A</i> subsection: Add <i>Buffer B</i>: 0.067 M dibasic potassium phosphate AND Line 1 of <i>Mobile phase</i>: Change Methanol and <i>Buffer</i> to: Methanol and <i>Buffer A</i>

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
							AND Line 4 of <i>Sample stock solution</i>: Change with the aid of 330 mL of <i>Buffer</i>, to a 1000-mL volumetric flask containing 50 mL of <i>Buffer</i>. to: with the aid of 330 mL of <i>Buffer B</i>, to a 1000-mL volumetric flask containing 50 mL of <i>Buffer B</i>.
PENTAZOCINE <i>Chemical Information</i>	USP36–NF31	4734	31-May-2013	1-Jun-2013	USP37–NF32	USP37–NF32	Line 1: Remove all chemical information.
POTASSIUM CHLORIDE EXTENDED-RELEASE CAPSULES	USP36–NF31	4838	31-May-2013	1-Jun-2013	USP37–NF32	USP37–NF32	Line 2: Change <i>Potassium stock solution</i> and <i>Standard preparations</i> —to: <i>Standard stock</i>

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
							<i>solution and Standard solutions—</i> AND Line 1 of <i>Procedure:</i> Change for <i>Procedure</i> in the Assay under <i>Potassium Chloride Oral Solution</i>. to: for <i>Instrumental conditions</i> and <i>Analysis</i> in the Assay under <i>Potassium Chloride Oral Solution</i>, except use Assay <i>preparation</i> instead of <i>Sample solution</i>.
SELEGILINE H <i>Dissolution</i> YDROCHLORI <711> DE TABLETS	USP36–NF31	5120	31-May-2013	1-Jun-2013	USP37–NF32	USP37–NF32	Line 3 of <i>Chromatographi c system:</i> Change Chromatograph the <i>Standard</i>

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
							<i>solution</i>, and record the peak responses. to: The flow rate is 1.0 mL/min. Chromatograph the <i>Standard solution</i>, and record the peak responses.
BIOTECHNOLOGICAL DERIVED AMINO ACIDS—AMINO ACID ANALYSIS NO ACID ANALYSIS ANALYSIS GENERAL PRINCIPLES	USP36–NF31	619	31-May-2013	1-Jun-2013	USP37–NF32	USP37–NF32	Change the section title <i>Method 6—Postcolumn DABS-Cl Derivatization General Principle</i> to: <i>Method 6—Precolumn DABS-Cl Derivatization General Principle</i>
VERAPAMIL HYDROCHLORIDE ORAL SUSPENSION Assay	USP36–NF31	5558	31-May-2013	1-Jun-2013	USP37–NF32	USP37–NF32	Line 2 of <i>Mobile phase</i>: Change 0.01 M to: 0.01 N

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
AMITRIPTYLIN IDENTIFICATION HYDROCHLORIDE TABLETS	USP36–NF31	2464	31-May-2013	1-Jun-2013	USP37–NF32	USP37–NF32	AND Line 7 of Assay preparation: Change 10-mL to: 100-mL Line 2: Change <i>Sample solution</i>: Nominally 0.01 mg/mL of amitriptyline hydrochloride in methanol from a suitable amount of finely powdered Tablets. Filter a portion of the solution, and use the filtrate for analysis. to: <i>Sample stock solution</i>: Nominally 0.1 mg/mL of amitriptyline hydrochloride in methanol from a

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
MINOCYCLINE Assay FOR INJECTION	USP36–NF31	4375	31-May-2013	1-Jun-2013	USP37–NF32	USP37–NF32	suitable amount of finely powdered Tablets. Filter a portion of the solution, and use the filtrate. <i>Sample solution:</i> Nominally 0.01 mg/mL of amitriptyline hydrochloride from <i>Sample stock solution</i> in methanol Line 2: Change <i>Mobile phase, Standard preparation, Resolution solution, and Chromatographic system</i> —Proceed as directed in the Assay under <i>Minocycline Hydrochloride</i>. to:

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
							<i>Mobile phase</i> —Prepare a mixture of 0.2 M ammonium oxalate, 0.01 M edetate disodium, dimethylformamide, and tetrahydrofuran (600:180:120:80). Adjust with ammonium hydroxide to a pH of 7.2, and pass through a filter of 0.5-μm or finer porosity. Make adjustments if necessary (see <i>System Suitability</i> under <i>Chromatography</i> <621>). <i>Standard preparation</i> —Dissolve an accurately weighed

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
							quantity of USP Minocycline Hydrochloride RS in water to obtain a solution having a known concentration of about 500 µg of minocycline (C₂₃H₂₇N₃O₇) per mL. Use this solution within 3 hours. <i>Resolution solution</i> —Transfer 10 mg of USP Minocycline Hydrochloride RS to a 25-mL volumetric flask, add 20 mL of 0.2 M ammonium oxalate, and swirl to dissolve. Heat on a water bath at 60° for 180 minutes, and

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
							allow to cool. Dilute with water to volume, and mix. <i>Chromatographic system (see Chromatography <621>)</i>—The liquid chromatograph is equipped with a 280-nm detector and a 4.6-mm × 25-cm column that contains 5-μm packing L1, and is maintained at a constant temperature of about 40°. The flow rate is about 1.5 mL per minute. Chromatograph the <i>Standard preparation</i>, and record the peak responses

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
							as directed for <i>Procedure</i>: the capacity factor, k', is not less than 5.0 and not more than 11.5; the tailing factor for the analyte peak is not less than 0.9 and not more than 2.0; and the relative standard deviation for replicate injections is not more than 2.0%. Chromatograph the <i>Resolution solution</i>, and record the peak responses as directed for <i>Procedure</i>: the relative retention times are about 0.7 for epiminocycline

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
							and 1.0 for minocycline; and the resolution, <i>R</i>, between epiminocycline and minocycline is not less than 4.6. AND Line 1 of <i>Procedure</i>: Change Proceed as directed for <i>Procedure</i> in the Assay under <i>Minocycline Hydrochloride</i>. to: Separately inject equal volumes (about 20 µL) of the <i>Standard preparation</i> and the Assay <i>preparation</i> into the chromatograph, record the

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Errata Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
OXYCODONE AND ACETAMINOPHEN TABLETS	IDENTIFICATION N/A. <i>Thin-Layer Chromatography</i>	USP36–NF31	4645	31-May-2013		1-Jun-2013	USP37–NF32	USP37–NF32	chromatograms, and measure the responses for the major peaks. Line 2 of <i>Sample solution</i> : Change in a mixture of methanol and water (4:1). to: in a 5-mL mixture of methanol and water (4:1).
POTASSIUM BICARBONATE AND POTASSIUM CHLORIDE EFFERVESCENT TABLETS FOR ORAL SOLUTION	<i>Assay for potassium</i>	USP36–NF31	4834	31-May-2013		1-Jun-2013	USP37–NF32	USP37–NF32	Line 2: Change <i>Potassium stock solution</i> and <i>Standard preparations</i> —to: <i>Standard stock solution</i> and <i>Standard solutions</i> —AND Line 1 of <i>Procedure</i> :

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
POTASSIUM CHLORIDE EXTENDED-RELEASE TABLETS	Assay	USP36–NF31 4841	31-May-2013	1-Jun-2013	USP37–NF32	USP37–NF32	Change for <i>Procedure</i> in the Assay under <i>Potassium Chloride Oral Solution</i>. to: for <i>Instrumental conditions</i> and <i>Analysis</i> in the Assay under <i>Potassium Chloride Oral Solution</i>, except use Assay <i>preparation</i> instead of <i>Sample solution</i>. Line 5: Change <i>Potassium stock solution</i> and <i>Standard preparations</i>— to: <i>Standard stock solution</i> and <i>Standard solutions</i>— Line 1 of

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
									<i>Procedure:</i> Change in the Assay under <i>Potassium Chloride Oral Solution</i> . to: for <i>Instrumental conditions</i> and <i>Analysis</i> in the Assay under <i>Potassium Chloride Oral Solution</i> , except use Assay <i>preparation 1</i> or Assay <i>preparation 2</i> instead of <i>Sample solution</i> .
ISOTRETINOIN CAPSULES	PERFORMANCE TESTS/ <711>/Test 4	Revision Bulletin (Official October 01, 2012)	Online	31-May-2013		1-Jun-2013	USP37–NF32	USP37–NF32	Line 2 of <i>Medium</i> : Change 4.5% (v/v) of Milloxid L (lauryl dimethyl amine oxide) to: 4.5% (v/v) of

Monograph Title Section		Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
TROLAMINE SALICYLATE	Assay	USP36–NF31	5499	31-May-2013	1-Jun-2013	USP37–NF32	USP37–NF32	lauryl dimethyl amine oxide
								Line 3 of <i>Chromatographic system:</i> Change L1 to: L7
NITRIC ACID	ASSAY/ <i>Procedure</i>	USP36–NF31	2107	31-May-2013	1-Jun-2013	USP37–NF32	USP37–NF32	Line 1 of <i>Sample solution:</i> Change To 2 mL of Nitric Acid in a tared, glass-stoppered conical flask add 25 mL of water. to: Weigh 2 mL of Nitric Acid in a glass-stoppered conical flask, and add 25 mL of water.
FOSPHENYTOIN SODIUM INJECTION	Assay	USP36–NF31	3680	31-May-2013	1-Jun-2013	USP37–NF32	USP37–NF32	Line 1 of <i>Assay preparation:</i> Change Transfer an

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
							accurately measured volume of the Injection, equivalent to about 300 mg of fosphenytoin, to: Transfer an accurately measured volume of the Injection, equivalent to about 300 mg of fosphenytoin sodium,
MOXIFLOXACIN OPHTHALMIC SOLUTION	USP36–NF31	4414	31-May-2013	1-Jun-2013	USP37–NF32	USP37–NF32	Line 4 of <i>Resolution solution</i> : Change 0.1 mg per mg and 0.001 mg per mg, to: 0.1 mg per mL and 0.001 mg per mL,
POTASSIUM BICARBONATE EFFERVESCE	USP36–NF31	4833	31-May-2013	1-Jun-2013	USP37–NF32	USP37–NF32	Line 2: Change <i>Potassium stock</i>

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
NT TABLETS FOR ORAL SOLUTION									<i>solution and Standard preparations—to: Standard stock solution and Standard solutions—AND Line 1 of Procedure: Change for Procedure in the Assay under Potassium Chloride Oral Solution. to: for Instrumental conditions and Analysis in the Assay under Potassium Chloride Oral Solution, except use Assay preparation instead of Sample solution.</i>
POTASSIUM	Assay	USP36–NF31	4840	31-May-2013		1-Jun-2013	USP37–NF32	USP37–NF32	Line 2: Change

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
CHLORIDE FOR ORAL SOLUTION							<i>Potassium stock solution and Standard preparations—</i> to: <i>Standard stock solution and Standard solutions—</i> AND Line 1 of <i>Procedure:</i> Change for <i>Procedure</i> in the Assay under <i>Potassium Chloride Oral Solution</i> . to: for <i>Instrumental conditions</i> and <i>Analysis</i> in the Assay under <i>Potassium Chloride Oral Solution</i> , except use Assay <i>preparation 1</i> or Assay <i>preparation 2</i>

Monograph Title	Section	Source Publication	Page Number	Errata Post Date	Errata Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
TACROLIMUS CAPSULES	ADDITIONAL REQUIREMENT S/USP Reference Standards <11>	USP36–NF31	5257	31-May-2013		1-Jun-2013	USP37–NF32	USP37–NF32	instead of <i>Sample solution</i>. Line 11 of USP Tacrolimus System Suitability Mixture RS: Change and tacrolimus 8-propyl analog (3S,4R,5S,8S,9E,12S,14S,15R,16S,18R,19R,26aS)-5,6,8,11,12,13,14,15,16,17,18,19,24,25,26,26a-hexadecahydro-5,19-dihydroxy-3-[(E)-2-[(1R,3R,4R)-4-hydroxy-3-methoxycyclohexyl]-1-methylvinyl]-14,16-dimethoxy-4,10,12,18-tetramethyl-15,19-epoxy-8-propyl-3H

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
							-pyri do[2,1-c][1,4]oxaazacycl otricosine-1,7,2 0,21-(4<i>H</i>,23<i>H</i>)-tetrone. to: and tacrolimus 8-propyl analog (3<i>S</i>,4<i>R</i>,5<i>S</i>,8<i>R</i> ,9<i>E</i>,12<i>S</i>,14<i>S</i> ,15<i>R</i>,16<i>S</i>,18<i>R</i> ,19<i>R</i>,26<i>aS</i>)-5,6,8,11,12,13 ,14,15,16,17,18, 19,24,25,26,26<i>a</i> -hexadecahydro -5,19-dihydroxy- 3-<i>{(E</i>)-2-<i>[(1R,3R,4</i> <i>R</i>)-4-hydroxy-3-m ethoxycyclohex yl]-1-methylvinyl }-14,16-dimetho xy-4,10,12,18-te tramethyl-15,19 -epoxy-8-pro- pyl-3<i>H</i> -pyri do[2,1-c

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
BIOTECHNOLOGY-DERIVED ARTICLES—AMINO ACID ANALYSIS	USP36–NF31	619	31-May-2013	1-Jun-2013	USP37–NF32	USP37–NF32	][1,4]oxaazacyclotricosine-1,7,20,21-(4<i>H</i>,23<i>H</i>)-tetrone. Change the section title <i>Method 6—Postcolumn DABS-Cl Derivatization</i> to: <i>Method 6—Precolumn DABS-Cl Derivatization</i>
VITAMIN E SPECIFIC TESTS/ <i>Acidity</i>	USP36–NF31	5579	31-May-2013	1-Jun-2013	USP37–NF32	USP37–NF32	Line 1 of <i>Sample</i>: Change 40 mg to: 1.0 g
FOSPHENYTOIN SODIUM	USP Reference standards <11>	USP36–NF31 3679	31-May-2013	1-Jun-2013	USP37–NF32	USP37–NF32	Line 6: Change C₁₄H₁₅NO₂ to: C₁₄H₁₃NO₂
MINOCYCLINE HYDROCHLORIDE ORAL SUSPENSION	Assay	USP36–NF31 4376	31-May-2013	1-Jun-2013	USP37–NF32	USP37–NF32	Line 2: Change <i>Mobile phase</i> and <i>Chromatographic system</i>

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
							—Proceed as directed in the Assay under <i>Minocycline Hydrochloride</i>. to: <i>Mobile phase</i>—Prepare a mixture of 0.2 M ammonium oxalate, 0.01 M edetate disodium, dimethylformamide, and tetrahydrofuran (600:180:120:80). Adjust with ammonium hydroxide to a pH of 7.2, and pass through a filter of 0.5-μm or finer pore size. Make adjustments if necessary (see <i>System Suitability</i> under <i>Chromatography</i> <621>).

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
							<i>Chromatographic system (see Chromatography <621>)</i>—The liquid chromatograph is equipped with a 280-nm detector and a 4.6-mm × 25-cm column that contains 5-μm packing L1, and is maintained at a constant temperature of about 40°. The flow rate is about 1.5 mL per minute. Chromatograph the <i>Standard preparation</i>, and record the peak responses as directed for <i>Procedure</i>: the capacity factor, k', is not less than 5.0 and not

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
							more than 11.5; the tailing factor for the analyte peak is not less than 0.9 and not more than 2.0; and the relative standard deviation for replicate injections is not more than 2.0%. Chromatograph the <i>Resolution</i> solution, and record the peak responses as directed for <i>Procedure</i>: the relative retention times are about 0.7 for epiminocycline and 1.0 for minocycline; and the resolution, <i>R</i>, between epiminocycline

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
OXYCODONE ASSAY/	USP36–NF31	4648	31-May-2013	1-Jun-2013	USP37–NF32	USP37–NF32	and minocycline is not less than 4.6. AND Line 1 of <i>Procedure</i>: Change Proceed as directed for <i>Procedure</i> in the Assay under <i>Minocycline Hydrochloride</i>. to: Separately inject equal volumes (about 20 µL) of the <i>Standard preparation</i> and the Assay <i>preparation</i> into the chromatograph, record the chromatograms, and measure the responses for the major peaks. Line 7 of

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
TEREPHTHAL <i>Procedure</i> ATE							<i>Analysis:</i> Change R_U = internal standard ratio (peak response of oxycodone/peak response of ethylparaben) from the <i>Standard solution</i> R_S = internal standard ratio (peak response of oxycodone/peak response of ethylparaben) from the <i>Sample solution</i> to: R_U = peak response ratio of oxycodone to ethylparaben from the <i>Sample solution</i> R_S = peak response ratio of oxycodone to ethylparaben from the

Monograph Title Section	Source Publication	Page Number	Errata Post Date Sort ascending	Errata Official Date	Target Errata Print Publication	Target Online Fix Publication	Description
							<i>Standard solution</i>

Pagination

- [First page](#) [« First](#)
- [Previous page](#) [< Previous](#)
- ...
- [Page 32](#)
- [Page 33](#)
- [Page 34](#)
- [Page 35](#)
- [Page 36](#)
- [Page 37](#)
- [Page 38](#)
- [Page 39](#)
- [Page 40](#)
- [Next page](#) [Next >](#)
- [Last page](#) [Last »](#)